

# Nothing Beats a Quickie™

The Quickie Q7 features 7000 Series Aluminium and ShapeLoc Technology that together shift the axis of the mobility world. Q7 offers great performance and sleek style at an aluminium price.



## Unparalleled Performance

- Design
  - Q7's minimalist ultra modern design optimises your style. Choose from a wide range of functional and contemporary options.
- Rigidity
  - Ovalised tubing in axle plate and caster housing decreases frame flex.
  - Increased frame bend strength results in better propulsion efficiency.
- Adjustability
  - Q7 keeps the full range of ergonomic adjustments and the strength of a rigid frame without welding. The result is an ultralight chair that can be fine tuned to your personal requirements.
- Lightweight Lifting / Transfer
  - Two ergonomic grab bars with rubber grip pads make both tasks easy and efficient.

## Unbelievably Lightweight

- Fully loaded configuration as low as 6.4 kg.
- 7000 Series Aluminium
  - Lightweight Aluminium that can be tempered for strength through ShapeLoc technology.



**ShapeLoc Technology is the process of converting an Aluminium frame into stronger, more rigid T6 Aluminium. The open frame is fabricated using 7000 series Aluminium then heat treated to convert the metal to T6 Aluminium.**

**The Advantage: Stronger metal = thinner tube walls = an unbelievably lightweight chair!**

## SPECIFICATIONS Quickie Q7

Overall Weight	8 kg (40 x 40. Std Configuration)
Weight Capacity	120.5 kg
Product Width	30.5 - 50.8 cm (12" - 20")
Product Length	30.5 - 50.8 cm (12" - 20")
Centre of Gravity	2.5 - 10.2 cm (1" - 4")
Back Height Adjustments	20.3 - 50.8 cm (8" - 20")
Front Seat Height	40.6 - 53.3 cm (16" - 21")
Rear Seat Height	33 - 50.8 cm (13" - 20")
Frame Angle	75°, 80°, 85°
Warranty	Limited/Lifetime

**QUICKIE Q7**

[www.QuickieQ7.com](http://www.QuickieQ7.com)

© Sunrise Medical  
114489 Rev. A



Sunrise Medical Pty Ltd  
Unit 10, 15 Carrington Road  
Castle Hill NSW 2154  
PO Box 8161  
BHBC NSW 2153  
Ph: (02) 8858 9191  
Fax: (02) 8850 4711  
[orders@sunrisemedical.com.au](mailto:orders@sunrisemedical.com.au)  
[www.sunrisemedical.com](http://www.sunrisemedical.com)

**QUICKIE®**

**Q7**



**Alternative to Titanium...  
Lighter Weight,  
Lighter Price**

**QUICKIE Q7**

