

Nothing Beats a Quickie®



folding family **QUICKIE® 2**

QUICKIE®
MANUAL CHAIRS

As the industry benchmark for adjustable folding chairs, the Quickie 2 series fits the needs of people through all stages of life. The Quickie 2 continues to be a market leader because of its innovative design, versatility, modularity, quality and ease-of-use for which Quickie is famous.

Quickie 2 Series Manual Folding Wheelchairs

The modular chair for all stages of life



Innovative Design

The Quickie 2 is the most popular ultra-lightweight folding frame wheelchair on the market, ideal for riders with an array of **interchangeable frame components** that can easily be adjusted as needs evolve. Our **True Fit program** provides a one-time free replacement (to the original owner within three years) for growth of the seat sling, back upholstery, footplates, hangers, back tubes, frames and cross braces. (Also provided for Quickie 2HP.)




High Performance Options

Add sporty features, titanium hardware that's ultra-strong and very light, and more...and you get the **Quickie 2HP** – the "High Performance" cousin of the Quickie 2. Designed for a more active user, it features slotted caster forks, plus a **wide selection of performance wheels** for a rigid chair feel in a tubular folding frame. A quick adjustment system allows for **rear seat height** and **camber adjustments** in one-third the time of conventional systems. A sliding seat rail adds even greater fine tuning to seat depth.

The Quickie 2 family continues to lead the adjustable, folding wheelchair segment based on design versatility, quality and ease-of-use. The customization options of the Quickie 2 help people of all ages **Get Out There!**

To learn more, contact your local dealer, call Sunrise Medical at **1-800-333-4000** or visit **www.sunrisemedical.com**.



folding family QUICKIE® 2

QUICKIE®
MANUAL CHAIRS

	O2 (K0005)	O2HP (K0005)
INNOVATIVE DESIGN		
Frame Style	Folding-Traditional	Folding-Traditional
Transport Weight	15.2 lbs	13.2 lbs
Standard Weight	27 lbs	25 lbs
Weight Capacity	250 lbs STD; 350 lbs HD	250 lbs STD
CUSTOMIZATION/ADJUSTABILITY		
Seat Width	11"- 22"	11"- 22"
Seat Depth	10"- 20"	10"- 20"
Back Angle Adjustment	DAB: max range: 78°- 120°	DAB: max range: 87°- 105°
Camber	0°- 7°	0°- 7°
Tension Adjustable Uph	Yes	Yes
Overall Width	24.5"	25.5"
CLINICAL APPLICATIONS		
Front Seat-To-Floor Height	17"- 22"	17"- 22"
Rear Seat-To-Floor Height	15.5"- 21"	15.5"- 21"
Footrest/Frame Angle	Swing Away	Swing Away
Back Type	Folding/Non Folding	Non Folding
Back Height	8.5"- 21"	8.5"- 20"
Center Of Gravity (COG)	.05"- 3.5"	.05"- 3.5"
LIFESTYLE OPTIONS		
Xtender-Power Assist	Yes	Yes
Front Suspension	Yes	Yes
Rear Suspension	No	No
Color Fusion	Yes	Yes
Warranty	Limited Lifetime	Limited Lifetime

