


Knee Scooter

Datasheet: 00146



The Knee Scooter - The new alternative to continued mobility after foot and ankle injuries or surgery!

Description

- The Knee Scooter is designed to offer mobility to those who have undergone ankle, foot or lower-limb surgery, or who have injured their leg. The unique design of the Knee Scooter allows the user to still easily move around, without hindering or jeopardising the healing process. Commonly recommended by Physio departments as the more functional alternative to crutches.

Features

- Adjustable height handles
- Adjustable height knee rest
- Dual locking brakes
- Quick adjustment for left or right leg
- 8" Wheels
- Front carry basket
- Aluminium powercoated steel frame

Applications

This Knee Scooter is designed for use:-

- After Foot, Ankle or lower Leg surgery
- In the case of Ankle/Foot strains or breaks

Specifications

Knee-rest Height	430 - 585 mm
Handle Height	790 - 1000 mm
Unit Weight	12 KG
MUW	130 KG



P: 5144 4955 **F:** 5144 6044
E: sales@accessrehabequip.com.au
W: www.accessrehabequip.com.au
156 York Street, Sale, Vic

'Providing Solutions for Life'