



# Celebrity<sup>®</sup> DX

3 and 4-wheel scooters



*Easily disassembles into 6 components*



*Infinite adjustable tiller*



*Easy Drive tiller with wraparound handles*

The Celebrity<sup>®</sup> DX delivers superb performance and ease of use with sleek, sporty styling for an overall package that impresses with elegant simplicity. Loaded with standard features like a deluxe seat with headrest and manual recline, and dual-bulb lower light for pathway illumination, the Celebrity DX is equipped to excite.



Tel: 03 9706 4611  
Fax: 03 9706 4622  
[www.pridemobility.com](http://www.pridemobility.com)

**Pride**<sup>®</sup>  
Mobility Products  
Australia Pty. Ltd.

## Specifications

Weight capacity:	159 kg (350 lbs.)
Maximum speed:	Up to 9.5 km/h (6 mph)
Ground clearance:	9 cm (3.5")
Turning radius:	3 wh: 104 cm (41.5") 4 wh: 150 cm (59")
Overall length:	3 wh: 117 cm (46") 4 wh: 119 cm (46.75")
Overall width:	61 cm (24")
<b>Tyres:</b>	
Front:	9 cm x 25 cm (3.5" x 10")
Rear:	10 cm x 25 cm (4" x 10")
Anti-tips:	Rear
Standard seating:	Deluxe seat with headrest and manual recline
Suspension:	3 wh: rear 4 wh: Front and rear
Battery charger:	Off-board, 5 amp
Per-charge range:	32 AH: Up to 45 km (28 miles) 40 AH: Up to 48 km (30 miles)
Batteries:	(2) 32 AH standard, (2) 40 AH optional
Battery weight:	11 kg (24.5 lbs.) each
Total weight w/o battery packs:	3 wh: 72 kg (159.5 lbs.) 4 wh: 77 kg (170 lbs.)
Heaviest piece when disassembled:	30 kg (65.5 lbs.)
Warranty	2 year limited

## Colour



Red

*"The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice".*

*"Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement".*

**Pride**  
Mobility Products  
Australia Pty. Ltd.

21 Healey Rd.  
Dandenong, Victoria 3175  
[www.pridemobility.com](http://www.pridemobility.com)

Tel: (03) 9706 4611  
Fax: (03) 9706 4622


AUS-SC-CELEBRITY-DX-3&4-06-02-06 ACN 088 609 661 ABN 66 088 609 661

## Benefits

- State-of-the-art rear suspension
- Top speed of up to 9.5 km/h (6 mph)
- Maximum weight capacity of 159 kg (350 lbs.)
- Ultra-simple, auto-latching lockup for quick disassembly/reassembly
- Easily accessible tie down points (for transport of unoccupied scooter)
- Front basket and rearview mirror standard
- Tight turning radius of 150 cm (59")
- Includes standard angle-adjustable dual bulb kerb light
- Front frame-mounted seat post for optimal stability
- High visibility reflectors
- Microprocessor-based controller
- Tiller rotation lock permits tiller to be used as a handle during transport

## Options

- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Walker holder



*Infinite adjustable tiller*



