

Shoprider

888SE

ASK FOR THE BRAND SHOPRIDER



SHOPRIDER[®]
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

888SE

SHOPRIDER[®]
www.shoprider.com.au *mobility*

The SHOPRIDER™ 888SE is packed with all the features of the larger models with the advantage of its compact size. The 888SE equipped with large batteries and has the following features, 260mm tyres for safety & comfort, fully rotating lockable seat, adjustable steering tiller and fold up arm rests. The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS

888SE

Length	1280mm
Width	600mm
Height	1070mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	15 Degrees
Maximum Load Single User*	136Kgs
Charger	4Amp
Motor Size	0.75Hp
Maximum Speed	8Kph
Battery Capacity (Minimum Amp)	2 x 50Amp
Maximum Range #	up to 40km
Wheel Size	260mm
Turning Radius	1500mm
Weight (w/o batteries)	76kg

FEATURES

High Back Seat
Adjustable Steering Tiller
Adjustable Width Arm Rests
Fold Arm Rests
Sliding Rotating Seats
Front Shopping Basket
Variable Seat height
Horn
Battery Level Warning Light
Front & Rear Bumper Pads
Front & Rear Suspension
Front Lights
Anti Tips Wheels
Comprehensive Warranty


Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING: Do not exceed 10kph on public roads or paths

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



SHOPRIDER[®]
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au