

Wheelie Walkers

Datasheet: 00142



Wheelie-Walker - 6" & 8" Wheels, Folding

Applications

These wheelie-walkers are commonly used in:-

- Nursing Homes
- Hospitals
- Aged Care Facilities
- Residential

Specifications

Name	Code	Seat Height	Wheel Size	Handle Height Range	Width	Weighs	Folding	MUW
Mini	M4051	480mm	6"	740-850mm	400mm	6.5 KG	Yes	100 KG
Medium	M4053	560mm	6"	780-890mm	610mm	6.5 KG	Yes	100 KG
Petite	M4052	510mm	8"	750-860mm	570mm	7.5 KG	Yes	100 KG
Medium	M4054	560mm	8"	810-920mm	610mm	7.5 KG	Yes	100 KG
Tall	M4055	610mm	8"	860-970mm	610mm	7.5 KG	Yes	100 KG
Push-Down	M4050	540mm	5"	910mm	610mm	5.5 KG	Yes	100 KG

Features

- Ergonomic, moulded hand-grips
- Dual action brakes
- Folding frame and pivoting backrest bar for ease of transportation
- Comfortable padded seat with collapsible bag underneath
- Lightweight - weighs 6.5 KG (6" Style) & 7.5 KG for (8")
- Available in Sparkling Red, Blue or Green
- Range of sizes to suit individual needs



P: 5144 4955 **F:** 5144 6044
E: sales@accessrehabequip.com.au
W: www.accessrehabequip.com.au
 156 York Street, Sale, Vic

'Providing Solutions for Life'