

Patient Trolley

Datasheet: 00183

Specifications

Code: M4102

Max Height	Min Height	Width:	Length:	MUW:
830 mm	37-51 mm	840 mm	2230 kg	250 kg



Applications

The Alrick 4000 MK II Patient Transfer Trolley is designed specifically as a day surgery bed and recovery trolley.

Description

A compact and manoeuvrable trolley designed to withstand a rigorous and demanding environment. Anti-Static Central Locking Castors take the strain off the carer by making the trolley extremely smooth and easy to manoeuvre, while the Trendelenburg feature allows the carer to regulate the patients' blood circulation during transfer between the operating room and ward.

Features

- Electric Hi/lo
- Electric Backrest
- Electric Forward & Reverse Trendelenburg
- Battery Backup
- Push Handle
- 150mm Anti-Static Castors
- Optional side rails and IV Pole



P: 5144 4955 **F:** 5144 6044
E: sales@accessrehabequip.com.au
W: www.accessrehabequip.com.au
156 York Street, Sale, Vic

'Providing Solutions for Life'