

Aquatic Wheelchair

Datasheet: 00184

Applications

This lightweight wheelchair is ideal for use at leisure centres, hydrotherapy pools and for transit in hospitals and nursing homes. Designed specifically for use in wet areas and can tolerate full submersion into water.

Description

Robust plastic frame enables this wheelchair to be used as a commode or transport chair. Fold away arms and flip up footrest assist with transfers.

Water Wheelchair constructed from high strength PVC piping and fittings.

Features

- Robust fixed plastic frame
- PVC coated polyester water proof fabric
- Self-propelled push rings on wheels for easy manoeuvrability
- User brakes on chair for greater control
- Rounded edge design on back and arm rest for greater comfort
- Can be weighted for children
- Frame rapidly empties and fills with water for stability
- Adjustable brakes
- Swing away arms for side transfers
- Folding footrest for front entry
- High quality stainless steel castors
- 600mm large rear wheels

Specifications

Code: M5180

Seat width	440mm, 500mm
Seat height	570mm
Chair height	990mm
Chair width	680mm
Chair length	1050mm
Weight	15kg
MUW	120kg
HD MUW	200kg



P: 5144 4955 **F:** 5144 6044
E: sales@accessrehabequip.com.au
W: www.accessrehabequip.com.au
156 York Street, Sale, Vic

'Providing Solutions for Life'